# Juvora

---

## Instructions for use


Device: JUVORA™ Dental Disc

an **Invivo** company

3	Dental Prosthesis Forms	44	Norsk
4	English	48	Polski
8	Dansk	52	Português
12	Deutsch	56	Română
16	ΕΛΛΗΝΙΚΑ	60	Slovenščina
20	Español	64	Suomi
24	Français	68	Svenska
28	Hrvatski	72	Türkçe
32	Italiano	76	Lietuvių k.
36	Magyar	80	Eesti
40	Nederlands	84	Latviski

1. Telescope

2. Attachment

3. Implant Supported  
Superstructures

### Intended use

---

For the manufacture of frameworks used in the production of:

- i) Full and partial removable dentures and overdentures;
- ii) Implant dentures and implant bars;
- iii) Anterior and posterior crowns, and posterior 3-unit bridges (maximum 1 pontic)

### Contraindications

---

Not to be used in bridges with 2 or more pontics.

### Precautions and warnings

---

Warning: Do not mill the frameworks too thin.

Follow the milling parameter guidance given in the JUVORA™ processing guide:

- The minimum thickness of secondary crowns and matrix housings should not be less than 0.7 mm
- The minimum thickness for “conventional” retention grids should not be less than 2 mm.

\*Note: Talk to the manufacturer of your machine about the processing of PEEK. Follow the instructions of your CAM and the manufacturer of your machine relating to the processing of JUVORA™ dental discs.

- Minimum size of horizontal retainer in lower jaw 2 mm x 8 mm
- Minimum size of horizontal band in upper jaw 2 mm x 8 mm
- In the event of any packaging damage or being unintentionally opened before use, do not use the product
- Do not repair or extend frameworks after milling or initial implantation.

### Material Information

---

JUVORA™ dental discs are manufactured from Invibio PEEK polymer

Form:	Solid
Odour:	Odourless.
pH (Value):	Not known.
Boiling Point (°C):	Not known.
Melting Point (°C):	343
Flash Point (°C):	Not known.
Auto Ignition Temperature (°C):	595
Oxidising Properties:	Not applicable.
Vapour Pressure(Pascal):	Not known.
Solubility (Water):	Insoluble.

## Processing Instructions

---

Detailed processing instructions are provided and must be observed by the end users (dental labs and technicians) when milling frameworks from the JUVORA™ dental discs. Dental labs that are supplied with the disc will be certified through Juvora Ltd, a process that involves assessment of their capability to follow the processing guidance to manufacture high quality frameworks and finished dentures where applicable using their existing equipment. This assessment will be completed by Juvora Ltd or their certified milling centres.

## Finishing

---

For finishing work, hard metal cross-cut burs are recommended. Plain cross-cut burs or special fine-toothed burs are particularly suitable. Further fine finishing of horizontal bands and retainers can be carried out with a silicon polisher. Before high-gloss polishing, Juvora Ltd recommends brushing on the polishing motor. The device can be polished to a high gloss with a goat hair or cotton wool buffing wheel. For the insides of telescope attachments, Juvora Ltd recommends finebrushes. Use fine (white) polishing paste.

## Production

---

JUVORA™ dental discs blanks can only be processed on suitable milling machines with appropriate milling programs. The milling head must be suitable for machining PEEK polymer. Juvora Ltd recommends suitable adaptation of the CAM software by the software manufacturer to allow milling of JUVORA™. For an exact fit, a 5-axis milling machine of stable construction is recommended.

Juvora Ltd also

recommends cooling the milling head for more accurate milling and to prolong its useful life.\*See foot note. One dental disc may be used to mill more than one framework.

**Note:** “Do not repair or extend the JUVORA framework and instead mill a new framework from the JUVORA dental disc

## Prosthesis removal from the dental disc

---

Remove the prosthesis carefully from the blank using a hard metal cutter and pay particular attention to the edges of the prosthesis.

## Safety Information

---

JUVORA Oyster White: **Warning!** Hazardous respirable dust may be formed when used. Do not breathe dust. Safety data sheet available on request. While machining the JUVORA™ dental discs the following safety precautions are recommended:

- Dust mask or dust extraction
- Personal protective equipment(eye protection, gloves)

## Incident reporting

---

In the event of a serious incident in relation to the device, the user and/or patient shall report the incident immediately to the manufacturer and their local competent authority.

## Storage information

---

The Juvora™ dental disc should be stored in dry storage and avoid exposure to direct sunlight. The PEEK polymer from which the Juvora™ dental disc is made is stable and can be stored for an extended period of time (10 year shelf life). It has a working temperature range from cryogenic up to 250°C and hence the storage temperature range for the Juvora™ dental disc is any ambient temperature and humidity.

## Additional Information

---

For additional information contact [info@invibio.com](mailto:info@invibio.com)

## Disposal considerations

---

### Waste treatment methods

Disposal should be in accordance with local, regional, state or national legislation.

### Additional Information

The European waste codes are recommendations based on the scheduled use of this product. For alternative uses and applications, other waste codes may be allocated under certain circumstances. 07 02 13- waste plastic, 07 02 99-waste not otherwise specified.

## Safety and Performance for JUVORA™ Dental Disc

The JUVORA™ Dental Disc is an alternative to metal for the manufacture of dental prostheses. From a safety perspective compared to metal, the JUVORA™ material does not induce metal allergies and there is evidence which shows there is a reduction in allergic reactions. The performance of the device shows a higher strength to weight ratio compared with metal prostheses. The device is light weight and displays lower stiffness being more closely aligned to that of natural bone and displays higher shock absorption (dampening effect).

## Clinical Benefits of Juvora™ Dental discs

Compared to metal prostheses, the JUVORA™ material shows an improvement in patient in-mouth comfort and quality of life

## Symbols in association with the label and ISO 15223

	Date of manufacture		Use by date
	Manufacturer		Batch code
	Authorised Representative		Keep away from sunlight
	Consult Instructions for use		Keep dry
	Medical device		Catalogue number

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