

COMMENTARY ON:

What is the Best Material for an Interbody Cage?

COMMENTARY AUTHOR: Michael Veldman – Invibio Biomaterial Solutions™

AUTHORS OF ORIGINAL PAPER: Heeren S. Makanji MD, Gregory D. Schroeder MD, Alexander R. Vaccaro MD PhD, Eve G. Hoffman MD

JOURNAL: Clinical Spine Surgery, 2020 May;33(4):137139. doi:10.1097/Bsd.0000000000000869

LEVEL OF EVIDENCE: Level 2a Peer-Reviewed, Systematic Review of Cohort Studies

Overview

What is the best material for an interbody cage? This relevant question was recently discussed by orthopedic surgeons at The Rothman Institute, Thomas Jefferson University Hospital, Phila, Pa., in an article published in *Clinical Spine Surgery*. Dr. Alexander Vaccaro and Dr. Eve Hoffman analyzed the literature for PEEK interbody devices for TLIF while Dr. Gregory Schroeder and Dr. Heeren Makanji provided their assessment of new 3D printed titanium cages.

Summary

The case for PEEK-OPTIMA™ Polymer

PEEK-OPTIMA polymer has been used for Interbody fusion for over two decades. It was chosen in the late 1990's for its material modulus, which closely matches cortical bone, biocompatibility, and radiolucency which allows for better fusion assessment on X-ray. The modulus of PEEK is close to cortical bone at 3.5 GPa while the modulus of solid titanium is 30 times higher at 100-110 GPa.¹ This means that a cage made from PEEK can have better stress distribution with bony endplates compared to a titanium cage. It also means that the graft material inside a PEEK cage sees more stress and should remodel into solid bone according to Wolff's law.^{1,2} The authors report that,

*"Stress shielding due to modulus mismatch is the proposed mechanism for the observed rate of subsidence of titanium anterior cervical cages of up to 26%."*³⁻⁵

The radiolucency of PEEK can provide more accurate fusion assessment on postoperative imaging or CT scan compared to titanium which is radiopaque.⁶ The primary disadvantage of PEEK cited by the authors is that it is chemically inert which results in less cell adhesion and surface bone integration.

Enhancing PEEK for bone on-growth

With the aim of addressing these perceived disadvantages of PEEK, recent advances in surface modifications for PEEK have emerged including:

- Titanium coatings
- Surface porosity
- Surface chemistry changes with the addition of hydroxyapatite

These surface enhancements have been shown to improve direct bone contact with PEEK devices.^{3,7} To address the perceived disadvantages of titanium, advances in 3D printing have allowed for more porous titanium devices with surface topography intended to provide bone ingrowth and lower cage stiffness to address material modulus.^{8,9} The authors conclude that,

"Unfortunately, prospective trials comparing 3D-printed porous titanium cages to PEEK cages are lacking in the literature. However, the current literature comprising biomechanical analysis, ovine models, and retrospective human studies support the continued use of PEEK cages in TLIF rather than titanium cages."

3D Printed Titanium Cage for TLIF

Traditional titanium cages have been used for osteointegrative properties, but the high modulus of elasticity and stiffness has led to subsidence and implant migration.¹⁰ To address the issue of subsidence, 3D printed porous titanium cages have been developed to decrease the stiffness of the cage and mitigate the probability of subsidence. Most of the published information on 3D printed cages has been on animal models, which have shown some promise for bone ingrowth.¹¹

However, there is a **lack of human clinical data comparing 3D printed cages to PEEK or Enhanced PEEK cages and the notch sensitivity of porous titanium has already resulted in one device recall.**

PEEK cages are superior to Titanium in transforaminal lumbar interbody devices

The authors looked at six published studies which clinically compared PEEK and titanium cages.^{4-6,12-14} Only one of the studies was prospective and looked at ACDF rather than TLIF.⁶ This prospective study showed improvement in clinical outcomes including Neck Disability Index, Oswestry Disability Index and Visual Analog Scale (VAS) in the support of PEEK at 7 year follow up. Only one of the two studies investigating TLIF cages collected data on patient reported outcomes and this study showed statistically significant decrease in VAS back pain scores at 12 months in the PEEK group compared with the titanium group.¹² Five out of the six studies were published prior to the availability of 3D printed titanium technology. There are no prospective studies comparing 3D porous titanium cages to PEEK cages.

The authors conclude that, until there is further data on new technologies, **the literature continues to overwhelmingly support the use of PEEK cages for TLIF.**

Commentary

The interbody fusion market has seen a number of new surface enhancement technologies over the past decade including 3D printed Titanium and other roughened or porous surface technologies. While promising conceptually, these technologies still have a clinical burden of proof to demonstrate similar fusion rates, clinical outcomes, low subsidence rates and imaging compatibility when compared to PEEK-OPTIMA polymer. ▲

ABOUT THE AUTHOR

Michael Veldman

Michael Veldman is the Global Strategic Marketing Manager for Invivio Biomaterial Solutions where he has worked since 2011. Prior to joining Invivio, he spent 2 years at Noalign Orthopedics where he led the development of a minimally invasive IM Nail and prior to that he spent 9 years in Product Development at Medtronic Spinal & Biologics. He received a Bachelor of Science degree in Mechanical Engineering from Christian Brothers University in 1998 and an MBA from the University of Memphis in 2002.



REFERENCES

1. Heary RF, et al. Elastic modulus in the selection of interbody implants. *J Spine Surg.* 2017;3(2):163-167.
2. Xiao Z, Wang L, Gong H, Zhu D. Biomechanical evaluation of three surgical scenarios of posterior lumbar interbody fusion by finite element analysis. *Biomed Eng Online.* 2012;11(1):31.
3. Rao PJ, Pelletier MH, Walsh WR, et al. Spine interbody implants: material selection and modification, functionalization and bioactivation of surfaces to improve osseointegration. *Orthop Surg.* 2014;6:81-89.
4. Chou YC, Chen DC, Hsieh WA, et al. Efficacy of anterior cervical fusion: Comparison of titanium cages, polyetheretherketone (PEEK) cages and autogenous bone grafts. *J Clin Neurosci.* 2008;15:1240-1245.
5. Niu CC, Liao JC, Chen WJ, et al. Outcomes of interbody fusion cages used in 1 and 2-levels anterior cervical discectomy and fusion: titanium cages versus polyetheretherketone (PEEK) cages. *J Spinal Disord Tech.* 2010;23:310-316.
6. Chen Y, Wang X, Lu X, et al. Comparison of titanium and polyetheretherketone (PEEK) cages in the surgical treatment of multilevel cervical spondylotic myelopathy: A prospective, randomized, control study with over 7-year follow-up. *Eur Spine J.* 2013;22: 1539-1546.
7. Walsh WR, et al. Does PEEK/HA Enhance Bone Formation Compared With PEEK in a Sheep Cervical Fusion Model? *Clin Orthop Relat Res.* 2016 DOI 10.1007/s11999-016-4994-x
8. Pelletier MH, Cordaro N, Punjabi VM, et al. PEEK versus Ti interbody fusion devices: Resultant fusion, bone apposition, initial and 26-week biomechanics. *Clin Spine Surg.* 2016;29:E208-E214.
9. McGilvray KC, Easley J, Seim HB, et al. Bony ingrowth potential of 3D-printed porous titanium alloy: a direct comparison of interbody cage materials in an in vivo ovine lumbar fusion model. *Spine J.* 2018;18:1250-1260.
10. Han CM, Lee EJ, Kim HE, et al. The electron beam deposition of titanium on polyetheretherketone (PEEK) and the resulting enhanced biological properties. *Biomaterials.* 2010;31:3465-3470.
11. Yang J, Cai H, Lv J, et al. In vivo study of self-stabilizing artificial vertebral body fabricated by electron beam melting. *Spine.* 2014;39: E486-E492.
12. Nemoto O, Asazuma T, Yato Y, et al. Comparison of fusion rates following transforaminal lumbar interbody fusion using polyetheretherketone cages or titanium cages with transpedicular instrumentation. *Eur Spine J.* 2014;23:2150-2155.
13. Cabraja M, Oezdemir S, Koepfen D, et al. Anterior cervical discectomy and fusion: comparison of titanium and polyetheretherketone cages. *BMC Musculoskelet Disord.* 2012;13:172-180.
14. Tanida S, Fujibayashi S, Otsuki B, et al. Vertebral endplate cyst as a predictor of nonunion after lumbar interbody fusion: comparison of titanium and polyetheretherketone cages. *Spine.* (Phila Pa 1976). 2016;41:E1216-E1222.



Invibio Ltd.

Victrex Technology Centre
Hillhouse International
Thornton-Cleveleys
Lancashire
FY5 4QD, UK
Tel: +44 (0) 1253 898 000
FAX: +44 (0) 1253 898 001

Invibio Inc.

300 Conshohocken State Road
West Conshohocken, PA
19428
USA
Toll Free: 866-INVIBIO (468-4246)
Tel: (484) 342-6004
Fax: (484) 342-6005

**For further information please email us at info@invibio.com
or visit our website at:**

► Invibio.com

Victrex plc and/or its group companies (Victrex plc) believes that the information contained in this document is an accurate description of the typical characteristics and/or uses of the product(s) and is based on information that we believe is reliable. However, it is provided for information only. It is not intended to amount to advice on which you should rely and should not be construed as, or used as a substitute for, professional medical advice or other professional or specialist advice. In particular, it is the customer's responsibility to thoroughly test the product in each specific application to determine its performance, efficacy, and safety for each end-use product, device or other application. Suggestions of product uses should not be taken as inducements to infringe any particular patent. Mention of a product in this document is not a guarantee of its availability.

Victrex plc reserves the right to modify products, specifications and/or packaging as part of a continuous program of product development. Victrex plc makes no warranties, express or implied, including, without limitation, a warranty of fitness for a particular purpose or of intellectual property non-infringement, including but not limited to patent non-infringement, which are expressly disclaimed, whether express or implied, in fact or by law.

Further, Victrex plc makes no warranty to your customers or agents and has not authorized anyone to make any representation or warranty other than as provided above. Victrex plc shall in no event be liable for any general, indirect, special, consequential, punitive, incidental or similar damages, or any damages for harm to business, lost profits or lost savings, even if Victrex has been advised of the possibility of such damages regardless of the form of action. The foregoing does not seek to affect any liability (including to individual consumers) which cannot be excluded or limited under any applicable law.

Supporting information is available on request for claims referenced in this document.

Victrex plc is the owner or the licensee of all intellectual property rights in the this document. All rights are protected by copyright and other intellectual property laws and treaties around the world. All rights reserved.

VICTREX™ and SHAPING FUTURE PERFORMANCE™ are trademarks of Victrex plc or its group companies.

INVIBIO™, JUVORA™, PEEK-OPTIMA™, BIOMATERIAL SOLUTIONS™ are trademarks of Victrex plc or its group companies.