

Manufacturing Trauma Fracture Fixation Implants: Partnering with Invibio

AUTHOR: Sherri Gambill – Invibio Biomaterial Solutions

Overview

Metal implants have been used for over 50 years with generally good outcomes, but in some fractures, complications related to non-unions, delayed unions, and implant failure continue to be a challenge, with overly stiff constructs as a reported risk factor.^{1,2} PEEK-OPTIMA™ Ultra-Reinforced is a composite material growing in popularity as an alternative to stainless steel and titanium for fracture fixation devices, and has seen clinical success in a variety of application areas.³⁻⁶

PEEK-OPTIMA Ultra-Reinforced combines the high performance material properties of PEEK-OPTIMA Polymer with the strength imparted by continuous carbon fibers. When manufactured into trauma devices, PEEK-OPTIMA Ultra-Reinforced enables semi-rigid fixation with improved fatigue and imaging properties over metal.⁷ Unlike metals, the carbon fiber PEEK offers the ability to tailor the mechanical properties of the implant without altering the geometry, offering increased design flexibility to meet the device requirements.


Partnering with Invibio enables medical device companies to expand their trauma product offerings beyond metal technologies at a fraction of the time and investment it would take to translate from a metal to composite solution internally. Medical device manufacturers can leverage Invibio’s expertise in composite technology, state-of-the-art tools, and dedicated manufacturing facility to bring new fracture management solutions to market.

Idea to Innovation – do more with fewer resources

Medical device manufacturers are experts in implant design, but converting those designs to composites requires time and money to build the knowledge, processing capabilities and supply chains internally. Invibio has invested heavily to build the capabilities required to deliver these components including: basic research, composite processing knowledge, and application-specific knowledge. A dedicated staff, assembled from the polymer/composite and medical device industries, works with partners from concept through production to overcome design and manufacturing challenges by offering comprehensive assistance, including design



Design for Manufacturing



Prototyping and Testing

for manufacturing, prototyping, testing, and regulatory support. Invibio’s state-of-the-art manufacturing facilities are operated under ISO 13485 certified quality management system and feature dedicated medical testing laboratories and a controlled manufacturing environment. These investments enable medical device companies to iterate quickly with low risk and decreased investment compared to developing on their own.

Components of any size and length can be produced, from small extremity plates to large distal femur plates. The process of plate creation is controlled by Invibio, from monomer through plate production. The process begins with the production of PEEK-OPTIMA Natural, which is then combined with carbon fibers into a tape. The tape is cut and compression molded into semi-finished components, and then finished to the customer’s design specifications and inspected.



Clinical Relevance – why change materials?

Locked plating is a significant advancement in fracture care resulting in improved patient outcomes for certain fractures.⁸ However, early reports of clinical success were followed by reports of clinical failures, which suggested that in some applications, the plate and screw construct may be too rigid, inhibiting the interfragmentary micromotion necessary to permit secondary healing by callus formation.^{9,10} Distal femur fractures are an often cited example where overly stiff locking plate constructs may lead to healing difficulties, with recent studies reporting non-union rates up to 20%^{2,11-15}

Strategies have been developed for reducing construct stiffness in three areas: (1) modifications in the surgical technique for existing plates and screws, (2) new screw designs, and (3) material advancements. Focusing on material advancements, studies have shown that a material with a lower elastic modulus may improve outcomes. More flexible titanium plates produce more callus and have fewer non-unions than stainless steel plates.^{2,9,14-16} PEEK-OPTIMA Ultra-Reinforced provides another alternative to the goal of more flexible fixation.

Testing demonstrates a plate produced from PEEK-OPTIMA Ultra-Reinforced can have reduced stiffness and greater fatigue strength than a titanium plate of the same geometry.⁷

Performance – Design Flexibility through Material

Carbon Fiber PEEK polymer plate stiffness and strength come not only from the plate geometry, but from the orientation of the carbon fibers throughout the plate, offering a huge array of choices to meet device specifications. This design flexibility is why carbon fiber devices have been adopted not only in medical devices, but in many advanced applications including the aerospace and automotive industries.

The plate geometry does not need to change in order to alter mechanical properties. In a 4-point bend test (per ASTM F382) of four identical generic distal femur plates, changes to the order of fiber orientation enabled a reduction in stiffness without a significant impact to the yield strength. In the example of Variant A to B, reducing stiffness by 12% resulted in a loss of yield strength of only 2% (ref. figure 1).⁷

Distal Femur Plate: 4-Point Bend Testing⁷




Figure 1

Conclusion

Invivo's investment in research and development, staff, and facilities to produce composite trauma plates enables medical device manufacturers to provide new options for treating traumatic injuries, with the potential for improved OR and point of care efficiencies, increased confidence to progress patients through recovery, and fewer and easier revisions.⁷

ABOUT THE AUTHOR

Sherri Gambill

Sherri (Wykosky) Gambill is currently Trauma Technology Manager at Invivo Biomaterial Solutions where she is responsible for product development. Previously, as Business Development Associate, she maintained relationships across client organizations as they adopted new biomaterials. Prior to Invivo, Sherri was a Product Development Engineer at DePuy Synthes and BD Ophthalmic Systems, where she designed and developed implants and instrumentation for orthopaedic trauma and glaucoma treatment. In 2006, Sherri received a Bachelor of Science (BS) degree in Bioengineering at the University of Pennsylvania in Philadelphia, PA, USA.



REFERENCES

- Henderson CE, Kuhl LL, Fitzpatrick DC, Marsh JL. Locking plates for distal femur fractures: is there a problem with fracture healing? *J Orthop Trauma*. 2011 Feb;25 Suppl 1:S8-S14. doi: 10.1097/BOT.0b013e3182070127.
- Rodriguez EK, Zurakowski D, Herder L, Hall A, Walley KC, Weaver MJ, Appleton Pt, Vrahas M. Mechanical construct characteristics predisposing to non-union after locked lateral plating of distal femur fractures. *J Orthop Trauma*. 2016 Aug;30(8):403-8. doi:10.1097/BOT.0000000000000593.
- Caforio M, Perugia D, Colombo M, Calori GM, Maniscalco P. Preliminary experience with Piccolo Composite™, a radiolucent distal fibula plate, in ankle fractures. *Injury*. 2014 Dec;45 Suppl 6:S36-S38. doi: 10.1016/j.injury.2014.10.020.
- Hak, DJ, Mauffrey C, Oliver T. Seeing is believing: treatment of proximal humerus fractures using a novel radiolucent implant and its effect on reduction accuracy, healing rate, and functional outcome. Paper #482. Presented at EFORT 2014 meeting.
- Guzzini M, Lanzetti RM, Luariello D, Morelli F, Princi G, Perugia D, Ferretti A. Comparison between carbon-peek plate and conventional stainless steel

plate in ankle fractures. A prospective study of two years follow-up. *Injury*. 2017 Jun;48(6):1249-1252.

- Rives AW. Distal Radius ORIF using a carbon-fiber reinforced PEEK volar palte for treatment of a left displaced intra-articular fracture with loss of height, radial inclination and volar tilt. Case Report. https://media.wix.com/ugd/cd37de_c6a5c6415aef4098a8ede8f476ddc9c0.pdf. Accessed May 18, 2017.
- Data on file at Invivo. Biomechanical performance is not indicative of clinical performance.
- Strauss EJ, Schwarzkopf R, Kummer F, Egol KA. The current status of locked plating: the good, the bad, and the ugly. *J Orthop Trauma*. 2008 Aug;22(7).
- Lujan TJ, Henderson CE, Madey SM, Fitzpatrick DC, Marsh JL, Bottlang M. Locked plating of distal femur fractures leads to inconsistent and asymmetric callus formation. *J Orthop Trauma*. 2010 Mar;24(3):156-62. doi:10.1097/BOT.0b013e3181be6720.
- Röderer G, Gebhard F, Duerselen L, Ignatius A, Claes L. Delayed bone healing following high tibial osteotomy related to increased implant stiffness in locked plating. *Injury*. 2014 Oct;45(10):1648-52. doi: 10.1016/j.injury.2014.04.018.
- Henderson CE, Kuhl LL, Fitzpatrick DC, Marsh JL. Locking plates for distal femur fractures: is there a problem with fracture healing? *J Orthop Trauma*.
- Ehlinger M, Dujardin F, Pidhorz L, Bonneville P, Pietu G, Vandenbussche E. Locked plating for internal fixation of the adult distal femur-influence of the type of construct and hardware on the clinical and radiological outcomes. *Orthop Traumatol Surg Res*. 2014 Sep; 100(5):549-54. doi: 10.1016/j.ost.2014.06.005.
- Ricci WM, Streubel PN, Morshed S, Collinge CA, Nork SE, Gardner MJ. Risk factors for failure of locked plate fixation of distal femur fractures: an analysis of 335 cases. *J Orthop Trauma* 2014 Feb;28(2):83-9. doi:10.1097/BOT.0b013e31829e6dd0.
- Henderson CE, Lujan TJ, Kuhl LL, Bottlang M, Fitzpatrick DC, Marsh JL. 2010 Mid-America Orthopaedic Association Physician in Training Award: healing complications re common after locked plating for distal femur fractures. *Clin Orthop Relat Res* 2011 Jun;469(6):1757-65. doi:10.1007/s11999-011-1870-6.
- Hoffmann MF, Jones CB, Sietsema DEL, Tornetta P 3rd, Koenig SJ. Clinical outcomes of locked plating of distal femoral fractures in a retrospective cohort. *J Orthop Surg Res*. 2013 Nov 27;8:43. doi:10.1186/1749-799X-8-43.
- Gaines RJ, Sanders R, Sagi HC, Haidukewych GJ. Titanium versus stainless steel locked plates for distal femur fractures. OTA abstract. 2008. Paper Number 55.

Copyright ©2017 Invivo Ltd. INVIBIO™, PEEK-OPTIMA™, INVIBIO BIOMATERIAL SOLUTIONS™ are trademarks of Victrex plc or its group companies. All rights reserved.



Invibio Ltd.

Victrex Technology Centre
Hillhouse International
Thornton-Cleveleys
Lancashire
FY5 4QD, UK

Tel: +44 (0) 1253 898 000
FAX: +44 (0) 1253 898 001

Invibio Inc.

300 Conshohocken State Road
West Conshohocken, PA
19428
USA

Toll Free: 866-INVIBIO (468-4246)
Tel: (484) 342-6004
Fax: (484) 342-6005

**For further information please email us at info@invibio.com
or visit our website at:**

► Invibio.com

Victrex plc and/or its group companies ("Victrex plc") believes that the information in this document is an accurate description of the typical characteristics and/or uses of the product or products, but it is the customer's responsibility to thoroughly test the product in each specific application to determine its performance, efficacy, and safety for each end-use product, device or other application. Suggestions of uses should not be taken as inducements to infringe any particular patent. The information and data contained herein are based on information we believe reliable. Mention of a product in this document is not a guarantee of availability.

Victrex plc reserves the right to modify products, specifications and/or packaging as part of a continuous program of product development. Victrex plc makes no warranties, express or implied, including, without limitation, a warranty of fitness for a particular purpose or of intellectual property non-infringement, including, but not limited to patent non-infringement, which are expressly disclaimed, whether express or implied, in fact or by law.

Further, Victrex plc makes no warranty to your customers or agents, and has not authorized anyone to make any representation or warranty other than as provided above. Victrex plc shall in no event be liable for any general, indirect, special, consequential, punitive, incidental or similar damages, including without limitation, damages for harm to business, lost profits or lost savings, even if Victrex has been advised of the possibility of such damages regardless of the form of action.

Supporting information is available on request for all claims referenced in this document.

Copyright ©2017 Invibio Ltd. INVIBIO™, JUVORA™ PEEK-OPTIMA™, INVIBIO BIOMATERIAL SOLUTIONS™ are trademarks of Victrex plc or its group companies. All rights reserved.