


# HYDROXYAPATITE (HA) OVERVIEW

 Physical and chemical structure similar to the mineral component of bone – apatite makes up 70% of bone mass




- OSTEOCONDUCTIVE bone chemically 'bonds' to it
- BIOCOMPATIBLE lack of local or systemic toxicity
- Most common and stable of the Calcium Phosphates used in bone grafting
- Used in dental applications since the 1970's and orthopedics since the 1980's where it is commonly used to coat titanium implants
- Highly crystalline forms are very stable and essentially non-resorbable
- Synthetic forms carry no risk of disease transmission

# BENEFITS OF HA IN THE PEEK-OPTIMA™ MATRIX

- HA available on ALL surfaces integrated, not coated
- Osteoconductive AND load-bearing
   HA alone is brittle with low fracture resistance
- No delamination the HA in PEEK-OPTIMA HA Enhanced is NOT used as a coating technology



Free HA Particles AT X6000 MAGNIFICATION



# EXAMPLES OF HA IN CLINICAL USE<sup>1</sup>

### **Bone Void Fillers**

- MASTERGRAFT® (Medtronic) (60% HA/40% βTCP)
- CERAMENT® (Bonesupport)
  (40% HA/60% α-calcium sulphate hemi-hydrate)
- Endobon® (Biomet) (100% natural HA)

### **HA Coated Dental Implants**

- Zimmer® MP-1® HA Technology, has been used clinically since 1996
  - Currently used on Tapered Screw-Vent® Implants

### **HA Coated Orthopedic Implants**

- ▶ DePuy Synthes' CORAIL® Hip System with HA coated stem
  - Over 30 years clinical history and 1.8 million patients worldwide<sup>2</sup>
  - 96.3% survivorship of the CORAIL stem at 25 years<sup>3</sup>

## HA Coated Bone Screws and Pedicle Screws

- Reform® HA Coated Pedicle Screw System (Precision Spine)
- OsteoTite Bone Screws for extremity fixation and HA coated screws of the Firebird®
   Spinal Fixation System (Orthofix)
- ▶ SI-LOK® Sacroiliac Fixation System (Globus Medical)



300 Conshohocken State Rd Invibio Ltd.

896-INVIBIO (468 4246) West Conshohocken, PA USA

FAX: +484 342 6005

Tel: +484 342 6004

Tel: +44 (0) 1253 898000

FAX: +44 (0) 1253 898001

United Kingdom

Invibio Ltd.

Lancashire FY5 4QD

Thornton Cleveleys,

Hillhouse International Technology Centre,

### For further information call us toll free at 866-INVIBIO or +44 (0)1253 898000 or please visit our website at:

### moo.oidivnl

*REFERENCES* 

- 1. Example products only. Invibio makes no claims or recommendations in relation to these products or their performance, clinical or otherwise.
- 2. CORAIL hip system clinical summary. DePuy Synthes. CA#DSEM/JRC/1015/0344(1) Issued:01/17
- 3. Vidalain JP. 25-year ARTRO Results: A Special Vintage from the Old World. The CORAIL Hip System: A Practical approach based on 25 years of experience. 2011;Chapter 4.2.134-101

be taken as inducements to infringe any particular patent. The information and data contained herein are based on information we believe reliable. Mention of a product in this document is not a guarantee of availability. customer's responsibility to thoroughly test the product in each specific application to determine its performance, efficacy, and safety for each end-use product, device or other application. Suggestions of uses should not Victrex plc and/or its group companies ("Victrex plc") believes that the information contained in this document is an accurate description of the typical characteristics and/or uses of the product or products, but it is the

limitation, a warranty of fitness for a particular purpose or of intellectual property non-infringement, including, but not limited to patent non-infringement, which are expressly disclaimed, whether express or implied, Victrex plc reserves the right to modify products, specifications and/or packaging as part of a continuous program of product development. Victrex plc makes no warranties, express or implied, including, without

Further, Victrex plc makes no warranty to your customers or agents, and has not authorized anyone to make any representation or warranty other than as provided above. Victrex plc shall in no event be liable for any

of such damages regardless of the form of action. general, indirect, special, consequential, punitive, incidental or similar damages, including without limitation, damages for harm to business, lost profits or lost savings, even if Victrex has been advised of the possibility

Supporting information is available on request for all claims referenced in this document.

Copyright © 2017 Invibio Ltd. INVIBIO™, PEEK-OPTIMA™, INVIBIO BIOMATERIAL SOLUTIONS" are trademarks of Victrex plc or its group companies. All rights reserved.