



PEEK POLYMER

FOR DRUG DELIVERY DEVICES

MICROPUMP WITH WEARABLE PATCH

Insulin delivery systems

WEARABLE PATCH – BENEFITS OF PEEK

- Biocompatible
- Sterilization resistant
- Artifact-free imaging for unobstructed views



MICROPUMP – BENEFITS OF PEEK

- High mechanical strength
- Miniaturization of parts
- Wear resistant
- Water resistant
- Telemetry and communication characteristics
- Stable electric insulation and thermal dissipation
- Compatible with a range of joining technologies

WHAT IS PEEK?

A high performance, non-resorbable thermoplastic polymer of the PAEK (Polyaryletherketone) family, available in a wide range of forms for use in different manufacturing processes.

WHY PEEK FOR DRUG DELIVERY APPLICATIONS?

- ✓ Miniaturized parts
- ✓ Connected devices
- ✓ Low-cost manufacturing
- ✓ Reliable and accurate dosing
- ✓ Patient-friendly for treatment adherence
- ✓ Ease of processing, injection molded or machined



METERED DOSE INHALER

For Asthma and Chronic Obstructive Pulmonary Disease (COPD)

BENEFITS OF PEEK

- Compressive strength and creep* resistance
- Good dimensional stability
- Resistant to chemicals
- No interaction with the API
- No additives or processing agents

OTHER APPLICATION EXAMPLES

- Electronic Drug Delivery Systems
- Injection pen / Auto-injector devices
- Primary packaging

*Supporting information available upon request

THE #1 PEEK EXPERTS

VICTREX™ AND ITS SUBSIDIARY INVIBIO BIOMATERIAL SOLUTIONS™



**FOR NON-IMPLANTABLE
DRUG DELIVERY APPLICATIONS**

VICTREX™ PEEK Polymer

- Integrated supply chain
- Security of supply
- Local technical support
- Millions of people rely on Victrex PEEK polymers everyday*



**FOR IMPLANTABLE
DRUG DELIVERY APPLICATIONS**

Invibio PEEK-OPTIMA™ Natural Polymer

- Comprehensive Master File
- Zero material-related recalls
- Extensive clinical evidence
- ISO-certified manufacturer
- Used in medical devices cleared globally:*
- US FDA 510(K), CE-marked and other regulatory bodies

#1 PEEK EXPERTS

40+ YEARS OF INNOVATION

- 500M+ Automotive gears
- 20K+ Aircraft
- 4B+ Mobile devices
- Ventilators

 www.victrex.com

20+ YEARS OF CLINICAL HISTORY

~9M PEEK OPTIMA DEVICES IMPLANTED WORLDWIDE

- Interbody Fusion Devices for spinal fusion
- Trauma Fracture Fixation plates
- Suture anchors for shoulder repair
- Dental Prosthetic implants

 www.invibio.com

Victrex plc and/or its group companies ("Victrex plc") believes that the information contained in this document is an accurate description of the typical characteristics and/or uses of the product or products, but it is the customer's responsibility to thoroughly test the product in each specific application to determine its performance, efficacy, and safety for each end-use product, device or other application. Suggestions of uses should not be taken as inducements to infringe any particular patent. The information and data contained herein are based on information we believe reliable. Mention of a product in this document is not a guarantee of availability.

Victrex plc reserves the right to modify products, specifications and/or packaging as part of a continuous program of product development. Victrex plc makes no warranties, express or implied, including, without limitation, a warranty of fitness for a particular purpose or of intellectual property non-infringement, including, but not limited to patent non-infringement, which are expressly disclaimed, whether express or implied, in fact or by law.

Further, Victrex plc makes no warranty to your customers or agents, and has not authorized anyone to make any representation or warranty other than as provided above. Victrex plc shall in no event be liable for any general, indirect, special, consequential, punitive, incidental or similar damages, including without limitation, damages for harm to business, lost profits or lost savings, even if Victrex has been advised of the possibility of such damages regardless of the form of action.

INVIBIO™, PEEK-OPTIMA™, INVIBIO BIOMATERIAL SOLUTIONS™ are trademarks of Victrex plc or its group companies.

VICTREX™ and the Triangle (Device), are trade marks of Victrex plc or its group companies.

DD-MS-INV-E-221-A (6/2020)

*Supporting information available upon request