


# JUVORA™ DENTAL DISC PRODUCT SHEET

JUVORA is manufactured from PEEK-OPTIMA™ and delivers a shock absorbing, high performance polymer alternative to zirconia or metal, for permanent, fixed and removable dental frameworks.

JUVORA™ can be manufactured with most conventional CAD/CAM systems:

- 4-5 axis milling systems that can accommodate 98mm stepped discs, or \*95mm slotted discs designed to accommodate Zirkonzahn® milling systems
- uses carbide milling tools
- can be milled using either dry or wet milling systems

Please refer to product IFU and Processing Guide for specific manufacturing and processing instructions

	Dimensions - 1 unit	SKU	
		Oyster White	Natural
	98mm x 16mm	540103	540113
	98mm x 18mm	540104	540114
	98mm x 20mm	540105	540115
	98mm x 22mm	540106	540116
	98mm x 25mm	540107	540117
	98mm x 30mm	540108	540118
	*95mm x 16mm	540109	540111
	*95mm x 20mm	540110	540112



The JUVORA™ Dental Disc is FDA and CE cleared. For more information please visit: [juvoradental.com/eifu](http://juvoradental.com/eifu)

For education and training please visit online  
<https://guide.juvoradental.com> or your app store and search for  
'JUVORA GUIDE'



Sales enquiries: [info@invibio.com](mailto:info@invibio.com)

The information contained herein is believed to be an accurate description of the typical characteristics and/or uses of Juvora product(s) and services. However, it is your ultimate responsibility to determine the performance, efficacy and safety of using Juvora product(s) and services for a specific application. Suggestions of uses should not be taken as inducements to infringe any particular patent or as a representation that the product is suitable for such uses. Juvora makes no warranties, express or implied, including without limitation, a warranty of fitness for a particular purpose or of intellectual property non-infringement, including, but not limited to patent non-infringement, which are expressly disclaimed, whether express or implied, in fact or by law. Furthermore, Juvora makes no warranty to your customers or agents, and has not authorized anyone to make any representation or warranty other than as provided above. Juvora shall in no event be liable for any general, indirect, special, consequential, punitive, incidental or similar damages, including without limitation, damages for harm to business, lost profits or lost savings, even if Juvora has been advised of the possibility of such damages, regardless to the form of action.